

## SMOOTH MUSCLE CELL MEDIUM

**Product Type:** Smooth Muscle Cell Medium (w/o FBS)  
**Catalog Number:** P60125 SF

国内流通品はFBSを含んでおりません。別途、FBSを10mL添加してご使用ください。

### **Product Description**

Smooth Muscle Cell Medium (SMCM) is a complete medium designed for optimal growth of normal human smooth muscle cells *in vitro*. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals and a low concentration of fetal bovine serum (2%). The medium is bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO<sub>2</sub>/95% air. The medium is formulated (quantitatively and qualitatively) to provide a defined and optimally balanced nutritional environment that selectively promotes proliferation and growth of normal human smooth muscle cells *in vitro*.

### **Components**

- 500 ml of Basal Medium
- 5 ml of Smooth Muscle Cell Growth Supplement (SMCGS).
- 5 ml of penicillin/streptomycin solution (P/S solution)

別途、10mLのFBSを添加してください。

### **Prepare for use**

Thaw SMCGS, FBS and P/S solution at 37°C. Gently tilt the SMCGS tube several times during thawing to help the contents dissolve. Make sure the contents of the supplement are completely dissolved into solution before adding to the medium. Rinse the bottle and tubes with 70% ethanol, and then wipe to remove excess. Remove the cap, being careful not to touch the interior threads with fingers. Add SMCGS, FBS and P/S solution into basal medium in a sterile field, mix well and then the reconstituted medium is ready for use. Since several components of this medium are light-labile, it is recommended that the medium not be exposed to light for lengthy periods of time. If the medium is warmed prior to use, do not exceed 37°C. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

### **Product Use**

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use *in vitro* diagnostic or clinical procedures.

株式会社ケー・エー・シー  
試験事業部  
〒661-0978 兵庫県尼崎市久々知西町2丁目1-20

(お問合せ先)  
TEL: 06-6435-9747 FAX: 06-6435-9748  
URL: <http://www.kacnet.co.jp/>  
E-mail: [shiyaku-info@kacnet.co.jp](mailto:shiyaku-info@kacnet.co.jp)