



REF: P60121

PERICYTE MEDIUM KIT

Product Type: Pericyte Medium kit

Catalog Number: P60121

国内流通品は FBS を含んでおりません。別途、FBS を 10mL 添加してご使用ください。

Product Description

Pericyte Medium (ECM) is a complete medium designed for optimal growth of normal human vascular pericytes in vitro. It is a sterile, liquid medium which contains essential and nonessential amino acids, vitamins, organic and inorganic compounds, hormones, arowth factors, trace minerals and a low concentration of fetal bovine serum (2%). The medium is bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO₂/95% air. The medium is formulated (quantitatively and qualitatively) to provide a defined and optimally balanced nutritional environment that selectively promotes proliferation and growth of normal human vascular pericytes in vitra

© Components

- 500 ml of Basal Medium
- 5 ml of Pericyte Growth Supplement (PGS).
- 5 ml of penicillin/streptomycin solution (P/S solution)

別途、FBS を 10mL 添加してください。

© Prepare for use

Thaw PGS, FBS and P/S solution at 37°C. Gently tilt the PGS tube several times during thawing to help the contents dissolve. Make sure the contents of the supplement are completely dissolved into solution before adding to the medium. Rinse the bottle and tubes with 70% ethanol, and then wipe to remove excess. Remove the cap, being careful not to touch the interior threads with fingers. Add PGS. FBS and P/S solution into basal medium in a sterile field, mix well and then the reconstituted medium is ready for use. Since several components of this medium are lightlabile, it is recommended that the medium not be exposed to light for lengthy periods of time. If the medium is warmed prior to use, do not exceed 37°C. When stored in the dark at 4°C. the reconstituted medium is stable for one month.

Product Use

THESE PRODUCTS ARE FOR RESEARCH USE ONLY. Not approved for human or veterinary use, for application to humans or animals, or for use in vitro diagnostic or clinical procedures.

株式会社ケー・エー・シー

試薬事業部

〒661-0978 兵庫県尼崎市久々知西町2丁目1-20

(お問合せ先)

TEL: 06-6435-9747 FAX: 06-6435-9748

URL: http://www.kacnet.co.jp/ E-mail: cs-info@kacnet.co.jpC