



3T3-L1 MEDIA COMPOSITIONS

3T3-L1 Adipocyte Medium (cat # F-AM-1-L1)	3T3-L1 Preadipocyte Medium (cat # F-PM-1-L1)
DMEM / Ham's F-12 medium (1:1, v/v) HEPES pH 7.4 Fetal Bovine Serum (FBS) Biotin Pantothenate Human insulin Dexamethasone Penicillin Streptomycin Amphotericin B	DMEM, high glucose HEPES pH 7.4 Bovine Calf Serum (BCS) Penicillin Streptomycin Amphotericin B
3T3-L1 Differentiation Medium (cat # F-DM-2-L1)	3T3-L1 Basal Medium (cat # F-BM-1-L1)
DMEM / Ham's F-12 medium (1:1, v/v) HEPES pH 7.4 Fetal Bovine Serum (FBS) Biotin Pantothenate Human insulin Dexamethasone Penicillin Streptomycin Amphotericin B Isobutylmethylxanthine PPAR γ agonist	DMEM/Ham's F-12 medium (1:1, v/v) HEPES pH 7.4 Biotin Pantothenate

NOTE:

All media except cat# PM-1-L1 contain 3.15g/L D-glucose.

PM-1-L1 contains 4.5g/L D-glucose.

All media are also available without serum and/or phenol red free.

Please inquire for custom media requests.

株式会社 ケー・エー・シー

試薬事業部

〒661-0978 兵庫県尼崎市久々知西町2丁目1-20

(お問合せ先)

TEL : 06-6435-9747 FAX : 06-6435-9748

URL : <http://www.kacnet.co.jp/>

E-mail : cs-info@kacnet.co.jp